

A general multi-criteria approach for iterative (re)optimization

Sebastiaan Breedveld, Pascal Storchi,
Marleen Keijzer and Ben Heijmen

RKF 2 February 2007

Introduction

Treatment planning is a *multi-criteria* problem.

Generally this problem is handled by manually adapting the importance factors for the organs involved, to meet treatment constraints.

Modern algorithms automatically optimize importance factors to meet the treatment constraints.

Introduction

Imposing constraints for such algorithms have the following problems:

- the constraints can be too *strict*
- the constraints can be too *loose*

Introduction

Imposing constraints for such algorithms have the following problems:

- the constraints can be too *strict*
- the constraints can be too *loose*

Resulting plan will be *suboptimal*, or even *unacceptable* if some constraints are not met.

Manually finding these constraints is very time-consuming.

We propose a method to automatically optimize treatment constraints.

Optimizing Constraints

Our method is to find a set of constraints with the following properties:


- the final plan satisfies all hard constraints
- the set is Pareto optimal

Optimizing Constraints

The process of finding constraints has been automatized.

A *constraint list* helps in sorting the importance of meeting a constraint.

A multi-criteria algorithm uses this list to optimize the constraints. This is *independent* of the underlying optimization algorithm (fluence map, beam angle, wedges, etc.)



Constraint list

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

constraint nr	volume	type	critical	objective	priority
1	PTV	DV	42.42 Gy	100%	0
2	PTV	Max		47.78 Gy	0
3	Body	Max		47.78 Gy	0
4	Bowel	DV	35 Gy	20%	1
5	Bladder	DV	40 Gy	40%	2
6	Colon	DV	40 Gy	20%	2
7	Bowel	DV	20 Gy	50%	3
8	Bladder	DV	20 Gy	75%	3
9	Colon	DV	20 Gy	30%	3
10	Body	DV	30 Gy	40%	4

Constraint list

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

constraint nr	volume	type	critical	objective	priority
1	PTV	DV	42.42 Gy	100%	0
2	PTV	Max		47.78 Gy	0
3	Body	Max		47.78 Gy	0
4	Bowel	DV	35 Gy	20%	1
5	Bladder	DV	40 Gy	40%	2
6	Colon	DV	40 Gy	20%	2
7	Bowel	DV	20 Gy	50%	3
8	Bladder	DV	20 Gy	75%	3
9	Colon	DV	20 Gy	30%	3
10	Body	DV	30 Gy	40%	4

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	0%	0
2	PTV	Max		47.78 Gy	0	0
3	Body	Max		47.78 Gy	0	0
4	Bowel	DV	35 Gy	20%	0%	1
5	Bladder	DV	40 Gy	10%	0%	2
6	Colon	DV	40 Gy	20%	0%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.12	0
3	Body	Max		47.78 Gy	44.79	0
4	Bowel	DV	35 Gy	20%	28.83%	1
5	Bladder	DV	40 Gy	10%	14.86%	2
6	Colon	DV	40 Gy	20%	8.30%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.12	0
3	Body	Max		47.78 Gy	44.79	0
4	Bowel	DV	35 Gy	20%	28.83%	1
5	Bladder	DV	40 Gy	15%	14.86%	2
6	Colon	DV	40 Gy	20%	8.30%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.12	0
3	Body	Max		47.78 Gy	44.79	0
4	Bowel	DV	35 Gy	20%	28.83%	1
5	Bladder	DV	40 Gy	15%	14.86%	2
6	Colon	DV	40 Gy	20%	8.30%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.15	0
3	Body	Max		47.78 Gy	44.83	0
4	Bowel	DV	35 Gy	20%	27.03%	1
5	Bladder	DV	40 Gy	15%	14.15%	2
6	Colon	DV	40 Gy	20%	8.26%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.15	0
3	Body	Max		47.78 Gy	44.83	0
4	Bowel	DV	35 Gy	28%	27.03%	1
5	Bladder	DV	40 Gy	15%	14.15%	2
6	Colon	DV	40 Gy	20%	8.26%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.15	0
3	Body	Max		47.78 Gy	44.83	0
4	Bowel	DV	35 Gy	28%	27.03%	1
5	Bladder	DV	40 Gy	15%	14.15%	2
6	Colon	DV	40 Gy	20%	8.26%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.15	0
3	Body	Max		47.78 Gy	44.83	0
4	Bowel	DV	35 Gy	27%	27.03%	1
5	Bladder	DV	40 Gy	15%	14.15%	2
6	Colon	DV	40 Gy	20%	8.26%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.23	0
3	Body	Max		47.78 Gy	45.08	0
4	Bowel	DV	35 Gy	27%	26.04%	1
5	Bladder	DV	40 Gy	15%	14.15%	2
6	Colon	DV	40 Gy	20%	8.21%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.23	0
3	Body	Max		47.78 Gy	45.08	0
4	Bowel	DV	35 Gy	26%	26.04%	1
5	Bladder	DV	40 Gy	15%	14.15%	2
6	Colon	DV	40 Gy	20%	8.21%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.76	0
3	Body	Max		47.78 Gy	46.32	0
4	Bowel	DV	35 Gy	26%	25.00%	1
5	Bladder	DV	40 Gy	15%	14.06%	2
6	Colon	DV	40 Gy	20%	8.11%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.76	0
3	Body	Max		47.78 Gy	46.32	0
4	Bowel	DV	35 Gy	24%	25.00%	1
5	Bladder	DV	40 Gy	15%	14.06%	2
6	Colon	DV	40 Gy	20%	8.11%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	48.03	0
3	Body	Max		47.78 Gy	46.84	0
4	Bowel	DV	35 Gy	24%	24.15%	1
5	Bladder	DV	40 Gy	15%	13.97%	2
6	Colon	DV	40 Gy	20%	8.09%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.76	0
3	Body	Max		47.78 Gy	46.32	0
4	Bowel	DV	35 Gy	25%	25.00%	0
5	Bladder	DV	40 Gy	15%	14.06%	2
6	Colon	DV	40 Gy	20%	8.11%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.76	0
3	Body	Max		47.78 Gy	46.32	0
4	Bowel	DV	35 Gy	25%	25.00%	0
5	Bladder	DV	40 Gy	14%	14.06%	2
6	Colon	DV	40 Gy	20%	8.11%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.78	0
3	Body	Max		47.78 Gy	46.50	0
4	Bowel	DV	35 Gy	25%	24.04%	0
5	Bladder	DV	40 Gy	14%	12.81%	2
6	Colon	DV	40 Gy	20%	8.06%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.78	0
3	Body	Max		47.78 Gy	46.50	0
4	Bowel	DV	35 Gy	25%	24.04%	0
5	Bladder	DV	40 Gy	14%	12.81%	2
6	Colon	DV	40 Gy	8%	8.06%	2

Example run

Typical constraint list for rectum cancer with a prescribed dose of 44.65 Gy.

nr	volume	type	critical	objective	met	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.78	0
3	Body	Max		47.78 Gy	46.50	0
4	Bowel	DV	35 Gy	25%	24.04%	0
5	Bladder	DV	40 Gy	14%	12.81%	2
6	Colon	DV	40 Gy	8%	8.06%	2


This process is repeated until none of the constraints can be fastened anymore.

Result for the rectum

constraint nr	volume	type	critical	objective	result obj.	priority
1	PTV	DV	42.42 Gy	100%	100%	0
2	PTV	Max		47.78 Gy	47.77 Gy	0
3	Body	Max		47.78 Gy	46.36 Gy	0
4	Bowel	DV	35 Gy	20%	24.3%	1
5	Bladder	DV	40 Gy	40%	10.0%	2
6	Colon	DV	40 Gy	20%	7.6%	2
7	Bowel	DV	20 Gy	50%	36.2%	3
8	Bladder	DV	20 Gy	75%	25.0%	3
9	Colon	DV	20 Gy	30%	12.6%	3
10	Body	DV	30 Gy	40%	17.8%	4

Optimization algorithm used is a fluence map optimization which iteratively adapts voxel-dependent importance factors in a quadratic objective function.

A result



Note: result is still subjective to the used optimization algorithm, but tight constraints enforce conformity.

A constraint list is very suitable for generalization of a treatment planning.

Can a single constraint list be used as a class solution?

constraint nr	volume	type	critical	objective	priority
1	PTV	DV	42.42 Gy	100%	0
2	PTV	Max		47.78 Gy	0
3	Body	Max		47.78 Gy	0
4	Bowel	DV	35 Gy	20%	1
5	Bladder	DV	40 Gy	40%	2
6	Colon	DV	40 Gy	20%	2
7	Bowel	DV	20 Gy	50%	3
8	Bladder	DV	20 Gy	75%	3
9	Colon	DV	20 Gy	30%	3
10	Body	DV	30 Gy	40%	4

Class solution: Rectum

- the constraint list was tuned on 4 patients
- this is verified on 4 other patients
- mean results are compared with the clinical plan (CadPlan)

		bowel	bladder	colon
tuning group	1	23%	18%	4%
	2	10%	5%	22%
	3	5%	11%	15%
	4	13%	47%	17%
verification group	5	11%	14%	5%
	6	9%	6%	10%
	7	10%	10%	-5%
	8	4%	18%	37%
average		10%	16%	13%
standard deviation		6%	14%	13%

Conclusion

- optimizing treatment constraints by a prioritized constraint list is a powerful tool for multi-criteria treatment planning
- a properly tuned constraint list makes it possible to achieve good results in an automated process